

NAR Simple Thoracostomy Kit

REV030322

Convenient solution for performing a simple (finger) thoracostomy...

The **NAR Simple Thoracostomy Kit** provides the essential equipment needed to perform a **simple (finger) thoracostomy**. Packaged in a compact and easy-to-open vacuum-sealed bag, the kit contains all necessary components to perform the procedure in one package. Contents include a 3-pack Chloraprep® swapstick for skin prep, a #10 PenBlade protected scalpel, sterile gauze, curved Rochester Pean hemostat to penetrate parietal pleura and dissect intercostal muscle, HyFin Compact Twin Pack Chest Seal for spontaneously breathing patients and a marker to write on chest seal as needed.



Simple Thoracostomy Kit contents:



Dimensions:

- Packaged: L 4.5 in. x W 7.5 in. x D 0.75 in.
- Weight: 5.28 oz

Simple Thoracostomy Kit

Item #	NSN#
83-0027	

Special Features:

- Compact kit that contains all necessary equipment to perform simple thoracostomy
- Includes 3-pack skin prep for performing the procedure bilaterally if needed
- HyFin Compact Vented Twin Pack Chest Seal to cover site if required

Educational Training Resources Available



Download
[NAR Simple Thoracostomy Kit Training Power Point Slides](#)



View online
[NAR Simple Thoracostomy Kit Product Training Video](#)



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®