

# RUF KIT

REV113018

***Slim profile trauma aid pack designed to carry a minimum of two complete sets of IFAK contents along with additional trauma supplies***

The RUF Kit from North American Rescue is a slim profile trauma aid pack designed to carry a minimum of two complete sets of IFAK contents along with additional trauma supplies. Designed with two tear away internal panels that hold standard IFAK contents and two tear away clear accessory pouches for additional supplies. The internal kit also has a shock cord panel designed to hold a compact, soft litter. The kit can be worn with a single bandolier style sling strap as well as comes with a waist strap that can be removed based on the user's needs. The RUF Bag is also available as a bag only (No gear or PVC Med Patch)

## Features:

- Slim profile trauma aid pack with two sets of IFAK components plus extra trauma supplies
- Can be worn with single bandolier style sling strap and with a detachable waist strap
- Two tear away internal panels designed to hold standard IFAK contents
- Two additional tear away clear accessory pouches for additional equipment
- Internal shock cord area to hold compact soft litter or additional equipment
- External MOLLE on sides of bag for external tourniquet pouches
- Large panel on exterior of kit for ID patches

## Specifications:

- Dimensions: H 18 in. x W 8 in. x D 2 in.
- Kit Weight: 6.3 lb
- Bag Weight: 2 lb

## RUF Kit

Level	Item#
Basic	80-0974
with Combat Gauze LE	80-0975
Advanced	80-0976
Bag only	80-0967W

BLK



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**

# RUF KIT

Can be worn with single bandolier style sling strap as well as with a detachable waist strap

REV113018



Every RUF KIT includes these core items:



## PLUS+ BASIC

also includes:




## COMBAT GAUZE LE

also includes:



## ADVANCED also includes:



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**

WITH A

**MISSION®**