

SAM IO FIELD STORAGE CASE

REV050223



The SAM IO Field Storage Case is a custom nylon carrying case that includes elastic loops for storing needles as well as attachment points for the driver once you've removed these items from their original packaging. The case features two compartments with a simple divider. The divider provides a mesh pouch on the under side for storing additional loose pieces. The lower compartment contains sufficient room for storing training items or additional supplies.

Features:

- Custom field storage case designed for the SAM IO
- Elastic loops for storing needles and the SAM IO Driver
- Two compartments separated by a divider with convenient mesh pouch

Dimensions:

- Case Closed: H 7 in. x W 11.25 in. x D 5.5 in.
- Weight: 3 lb

SAM IO FIELD STORAGE CASE

COLOR	ITEM #	NSN #
	30-0254	



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®