

TCCC CLS KIT SKILLS SET

REV051723

For initial TCCC CLS care training

The TCCC CLS Kit Skills Set provides five training CLS Kits for student use during trauma lanes or scenarios. Each CLS Kit contains the necessary equipment for initial TCCC – CLS care in the standard Army CLS bag with the standard Army CLS packing list. The set contains 5 pre-packed CLS kits in a hard, wheeled case that allows for easy transport and storage.

The CLS Skills Kit is engineered for student use during trauma lanes or scenarios. Each CLS Skills Kit contains the necessary equipment for initial CLS care in the standard Army CLS bag with the standard Army CLS packing list.

Features:

- Contains pre-packed 5 CLS kits for use by students during trauma lanes or training scenarios
- CLS kits (also available for separate purchase) packaged in a self-contained, wheeled, hard case for easy transport and storage

Hard Case Features

- Retractable extension handle
- In-line skate style wheels for easy transport
- Padlock holes include stainless steel padlock protector clips
- Cushion grip handles on three sides for comfortable lifting and carrying
- Waterproof and dust tight design
- Ambient pressure equalization valve
- Molded in hinge latch system
- Trigger release latch system
- Resistant to UV, solvents, corrosion and fungus
- Tan color

Specifications:

- External (Case): L 32.5 in. x W 21 in. x D 13 in.
- Internal (Case): L 29 in. x W 18 in. x D 10 in.
- Weight: 68 lbs



TCCC CLS Skills Kit also available for individual sale.



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

TCCC CLS KIT SKILLS SET

Option	Item#	NSN#
Skills Set	80-1212	----
Ind. Skills Kit	80-1199	----