

NAR QUIKLITTER™ *Non-Rigid Litter*

REV051923



Lowest cost solution for safe & effective transport

North American Rescue's QuikLitter™ is a cost effective emergency carry litter designed for rapid deployment for the movement of casualties and is the ideal disposable litter for mass casualty incidents, disasters, active shooter events, as well as for everyday patient movements.

Constructed of a water-resistant, woven polyethylene that is latex free and rated for up to 500 lb, the NAR QuikLitter™ features ten (10) carry handles with four (4) on each side and one (1) larger handle at either end of the litter to provide multiple carry points.

Features:

- Cost effective emergency carry litter
- Features ten large carry handles for ease of patient handling/transport
- Load capacity tested to 500 lb
- Latex free
- Water-resistant

Standard Size Specifications:

- Litter Bed: L 74 in. x W 28 in.
- Packaged: L 11 in. x W 7.75 in. x D 2 in.
- Weight: 13.3 oz
- Payload: 500 lb

NAR QuikLitter™



Size	Item#	NSN#
Standard	60-0061	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®