

# NAR QUIKLITTER™ XL *Non-Rigid Litter*

REV051923



## Lowest cost solution for safe & effective transport

North American Rescue's QuikLitter™ XL is a larger, cost-effective emergency carry litter designed for rapid deployment during casualty movement and is the ideal disposable litter for mass casualty incidents, disasters, active shooter events, as well as for everyday patient movements.

This multi-use, water-resistant, disposable evacuation platform is 12 in. wider and 6 in. longer than our standard QuikLitter™ and features 14 hand-holds (six on each side and one larger handle at each end) plus two built-in patient retention straps for easier and safer patient handling.

The QuikLitter™ XL is rated to 1,000 lb payload, offering more room, easier handling and overall safer patient transport.

### Features:

- Cost effective emergency carry litter
- Features (14) fourteen large carry handles with a load capacity tested to 1000 lb
- Includes (2) two built-in patient retention straps for easier and safer patient handling
- Latex free, water-resistant

### XL Size Specifications:

- Litter Bed: L 80 in. x W 40 in.
- Packaged: L 11 in. x W 7.75 in. x D 3 in.
- Weight: 2.2 lb
- Payload: 1,000 lb

NAR QuikLitter™ XL



Size	Item#	NSN#
XL	60-0065	



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**