

# NAR K-9 LITTER & Carrying Strap Accessories

REV092922



**Designed for the safe & efficient movement of an injured canine...**

The NAR K9 Litter is a lightweight, compact litter designed for the movement of injured canines. The vacuum sealed nylon litter can be packed to store in a cargo pocket or pouch and can be quickly deployed. The litter has leg cut outs to place the canine in the litter for safe movement. The K9 Litter also has six large and four small carry handles for balanced carrying and allows for one to four person carries. The resealable bag includes the NAR K-9 Litter and a K-9 Carry Strap & 2 Carabiners to aid in securing an over the shoulder carry of an injured canine with the NAR K-9 Litter. The maximum weight capacity for each K-9 litter is 700 lbs.

### Features:

- Lightweight, compact litter for the movement of an injured canine
- Can be carried in cargo pocket or pouch
- Leg cut outs to safely carry the canine in the litter
- Six large and four small carry handles that allow for one to four person carries
- Includes the K-9 Carry Strap & 2 Carabiners to aid in securing an over the shoulder carry of an injured canine using the NAR K-9 Litter

### Dimensions:

- Litter: W 35 in. x H 48 in.
- Adjustable Strap: L 9 to 54 in. x W 2 in.
- Carabiner open gate clearance: 27 mm
- Packaged: W 5.5 in. x H 10 in. x D 2.5 in.
- Weight: 1.38 lb
- Weight Capacity: 700 lbs

### NAR K-9 Litter & Accessories

| Options                                  | Item #   | NSN# |
|--|----------|------|
| K-9 Litter Only                          | 60-0074W |      |
| Carry Strap & Carabiners                 | ZZ-1004  |      |
| K-9 Litter with Carry Strap & Carabiners | 60-0076  |      |

### Carry Strap Accessory features:

- 2 in. wide adjustable seatbelt strap with durable plastic buckle to aid in securing an over the shoulder carry of an injured canine using the NAR K-9 Litter
- Includes 2 full-size carabiners with snag-resistant keylock nose design, large and ergonomic angled gate openings, machined frame ensures tight tolerances and smooth operation



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**

WITH A

**MISSION®**