

NAR RESCUE SLEEVE II™ & Accessories

REV022522



NOTE: Spine board sold separately.
Hoisting straps & carabiners not included.

BLK
ORG

SEL#: 09ME-04-SPIN
SEL#: 03WA-02-BORD

The Ultimate In Patient Immobilization...

The NAR Rescue Sleeve II™ is the industry's ultimate immobilization system supporting both vertical and horizontal hoisting. Simple and secure, the Rescue Sleeve II™ gives confidence to operators and comfort to patients. Used extensively in the field and refined for enhanced stability and overall functionality, it is standard issue to all U.S. Navy ships.

Features:

- Vertical and horizontal lifting-strap attachments
- Removable head-securing blocks
- Adjustable straps secure head and chin
- Hoist-cable storage space
- Padded head flaps
- Backboard compartment
- Padded carrying handles
- Easy to apply
- Replaceable parts are also available:
Head Block (#50-0042)
Head and Chin Straps (#50-0001 | NSN# 6515-01-477-8921)
Carrying Case (# 50-0002 | NSN# 4240-01-477-9408)
- Standard issue to all U.S. Navy ships



Specifications:

- Folded:
L 25.5 in. x W 13 in. x D 8 in.
- Open:
L 73 in. x W 18 in. x D 2 in.
- Payload: 1000 lb (453.6 kg)
- Weight: 16 lb

NAR Rescue Sleeve II

COLOR	ITEM#	NSN#
Black	50-0011	6530-01-531-6649
Orange	50-0012	6530-01-477-8525



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®