

# PERSONAL SAFETY CUTTERS *By Benchmade®*

REV012418

7 Hook



8 Hook

**BENCHMADE®**



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

## **Quick and Reliable During Emergency Rescues...**

Hook cutters designed to cut through a wide range of materials quickly and reliably during emergency rescues.

The larger handle of the 7 hook design works well with gloved hands and the coated handle allows for better grip. MOLLE compatible.

All of Benchmade's safety cutters are made along side their knives in the U.S.A

### **Features:**

- Saves time – more than 1.5 minutes faster than traditional trauma shears when cutting 2 articles of clothing.
- Extensively tested and widely used by the US military for medical and egress purposes
- Maintains sharpness for hundreds of cuts (up to 500+ cuts through seatbelt test material)
- Both hooks share the same 440C cutting surface and opening with varying handles to provide a more compact carry or more leverage
- Includes a MOLLE-compatible carrying sheath

### **8 Hook Specifications:**

- Open Length: 6.5 in.
- Weight: 1.5 oz

### **7 Hook Specifications:**

- Open Length: 4.3 in.
- Weight: 1.5 oz

### **Personal Safety Cutters By Benchmade®**

Model	Item#	NSN#
7 Hook	ZZ-0652	4240-01-543-9618
8 Hook	ZZ-0250	4240-01-568-3219

**PRODUCTS**  
WITH A  
**MISSION®**

# PERSONAL SAFETY CUTTERS *By Benchmade®*

REV012418



**MOLLE/PALS  
COMPATIBLE**  
*(hard sheath included  
with 7 hook only)*



**MAINTAINS SHARPNESS  
FOR HUNDREDS OF CUTS**  
*(up to 500+ cuts through  
seatbelt test material)*



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**