

# EMS RAPID DEPLOYMENT BAG *(Bag Only)*

REV010818



## **Highly visible, easily transportable and convenient...**

The EMS Rapid Deployment Bag is a bright orange bag with reflective trim for high visibility. Constructed of 1000 Denier Nylon, this rugged multiple load carriage configurable bag is easily transportable and highly identifiable.

The large central compartment and front trauma shear storage (trauma shears are sold separately) offer the user a versatile and mobile platform for transport.

### **Features:**

- Bag Only
- Sand & water resistant zipper with silent pulls
- Multiple-load carriage configuration capabilities: fanny pack, over the shoulder, hand-carried, attachments for MOLLE/PALS style webbing
- Constructed of 1000 Denier Nylon

### **Bag Specs:**

- H 7.75 in. x W 13 in. x D 9 in

### **EMS Rapid Deployment Bag**

Color	Item#	NSN#
ORG	80-0152T	



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**