

FUJI IX DENTAL CEMENT

REV022522

Component of the Dental Emergency Response Kit

Fuji IX GP is packable posterior glass ionomer restorative cement. This alternative to composites is perfect for geriatric and pediatric posterior restorations, final restorations, intermediate restorative, core material and long-term temporary restoration.

Specifications:

- Fluoride releasing
- Radiopaque
- Easy to handle and place
- Cures extremely hard and is very wear resistant
- Chemically bonds to tooth structure
- Tooth-like coefficient of thermal expansion

Dimensions:

- L 2.75 in. x W 2 in. x D 3 in.
- Weight: 0.25 lbs

Fuji IX Dental Cement

Rx

ITEM #	QUANTITY	NSN #
ZZ-1422	1	--



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®