

3-IN-1 DISPOSABLE ORAL EXAM PACK *Sterile*

REV022522

Component of the Dental Emergency Response Kit

This Intra-oral exam pack includes 1 disposable dental mirror with anti-fog coating and 1 disposable double ended probe and explorer. The probe is marked 3-6-9 MM universal style, tapered with round tip for readability and patient comfort. Explorer is made of high-grade stainless steel.

Specifications:

- 3-in-1 Intra-oral exam pack
- Sterile single use and disposable
- 1 dental mirror and 1 disposable double ended probe and explorer
- Probe marked with 3-6-9 MM universal style
- Explorer is made of high-grade stainless steel

Dimensions:

- L 9 in. x W 2.25 in.
- Weight: 0.03 lbs



3-in-1 Disposable Oral Exam Pack ☒ Rx **STERILE**

ITEM #	QUANTITY	NSN #
ZZ-1413	5 Needles	--
ZZ-1432	50 Needles	--



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®