

WHEELED LITTER CARRIER

REV052523

Lightest wheeled litter cart on the market for single patient transport...

The Wheeled Litter Carrier from North American Rescue is the lightest wheeled litter cart on the market for single person transport to & from MedEVAC vehicles, aircraft, etc. With its secure locking clamp system, this carrier can hold the litter firmly in place while allowing for different litter widths.

Compatible with Talon, Raven, MedEvac & T3 litters, the Wheeled Litter Carrier also performs as a litter stand in the field due to its unique wheel locking system. The carrier also is available with a rugged, stackable storage case with wheels for easy storage and transport.

The Wheeled Litter Carrier is quick and easy to deploy or disassemble for storage. High strength materials allow for a payload of 550 lb, and rugged, high-strength wheels are suitable for most environments and terrains.



NOTE: Litter sold separately



Carrier with Case is also available for sale ▶



Features:

- The lightest wheeled litter cart on the market
- Adjustable for different litter widths
- Compatible with Talon, Raven, MedEvac & T3 litters
- Also performs as litter stand in the field
- Easily foldable for transport & storage

Specifications:

- Wheeled Litter Carrier
Deployed: L 53.97 in. x W 27 in. x H 32.24 in.
Folded: L 41.73 in. x W 24.45 in. x H 9.37 in.
- Weight: 42 lb (19.05 kg)
- Payload: 550 lb (250 kg)
- Case Dimensions:
L 48 in. x W 29.5 in. x H 12.5 in.
Weight (without Carrier): 56 lb (25.4 kg)
Wheeled Litter Carrier with Case
Weight: 98 lb (44.4 kg)

NAR Wheeled Litter Carrier

Options	Item#	NSN#
Carrier only	60-0041	6530-01-663-2308*
Carrier with Case	60-0063	
Carrier with Case	60-0063-EU	

*NATO NSN



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®