

TOFFLEMIRE MATRIX BAND *Pack of 36*

REV022522

Component of the Dental Emergency Response Kit

The Tofflemire Matrix Bands are made with surgical grade stainless steel band is carefully stamped for consistently smooth, burr-free edges that keep your patient comfortable and reduce the risk of gingival trauma.

Specifications:

- Surgical grade stainless steel
- Smooth, bur free edges for patient comfort
- Single use and disposable
- Package of 36

Dimensions:

- L 6 in. x W 4 in.
- Weight: 0.01 oz



Tofflemire Matrix Band



ITEM #	QUANTITY	NSN #
ZZ-1417	36 Pack	--



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®