

TCCC IFAK SKILLS SET *All Combatant*

REV093021

Train up to 25 Students at Individual Skill Stations...

TCCC IFAK Skills Set – All Combatant: Provides the necessary equipment to perform required tasks when training at individual skill stations. Includes all the necessary equipment to instruct up to twenty-five (25) students in the treatment of the leading causes of potentially preventable death in the tactical environment. Contents are packaged in a self-contained, wheeled, rigid case with internal dividers to keep the equipment organized and easy to inventory. Addresses extremity hemorrhage, wound packing, chest seal application, needle chest decompression, basic airway adjuncts, splinting, hypothermia prevention & management, and eye injuries.

Features:

- The necessary training equipment for initial TCCC care using the IFAK Skills Set Kit- All Combatant
- Includes 25 Individual First Aid Kits (IFAK) each packaged in an easy to attach modular carrying bag that is compact and designed to be worn on the assault vest or body armor via MOLLE-style attachment points

Dimensions:

- Each IFAK: H 7 in. x W 4.5 in. x D 2.75 in.
Weight: 2 lb 8 oz
- Hard Storage Case:
External – L 32.5 in. x W 21 in. x D 13 in.
Internal – L 29 in. x W 18 in. x D 10 in.
- Total Weight: 70 lb

TCCC IFAK Skills Set – All Combatant

| Item# | NSN# |
|---------|------|
| 80-0969 | |



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
— WITH —
MISSION®

TCCC IFAK SKILLS SET *All Combatant*

REV093021



Each of the 25 TCCC IFAK Skills Kits contains:

- 1 x Nylon Bag - Coyote (Tactical Operator Response Kit)
- 1 x C-A-T® Tourniquet - Blue Trainer
- 1 x 6 in. ETD (Flat)
- 1 x NAR S-Rolled Gauze
- 1 x Combat Gauze Trainer
- 1 x 28F Naso with Lube
- 1 x HyFin® Trainer Twin Pack
- 1 x ARS® 14 ga x 3.25 in.
- 1 x Pair Bear Claw Nitrile Gloves - Lg
- 1 x Permanent Marker - Black, Lg
- 1 x Combat Casualty Card
- 1 x Roll 2 in. Surgical Tape



85-1958 TCCC IFAK Skills Kit – All Combatant is also available for individual purchase

Hard Case Specifications:

- Retractable extension handle
- In-line skate style wheels for easy transport
- Padlock holes include stainless steel padlock protector clips
- Cushion grip handles on three sides for comfortable lifting and carrying
- Waterproof and dust tight design
- Ambient pressure equalization valve
- Molded in hinge latch system
- Trigger release latch system
- Resistant to UV, solvents, corrosion and fungus
- Tan color



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A
MISSION®