

# NAR QUICK DETACH PANELS *3 size options*

REV103117



***Allows the operator the fastest access to their kit during the critical moments after wounding...***

The new line of NAR Quick Detach Panels from North American Rescue allows the operator the fastest access to their kit during the critical moments after wounding. Utilizing Mil-spec hook & loop and heavy-duty metal snaps to ensure that the detachable front panel remains secure in the most austere of conditions, the Quick Detach Panel can be used with any Mil-spec standard Molle/PALS system. Three versions of the panel designed to compliment all of NAR IFAK's are available including a NEW horizontal mount option

**Features:**

- **Small** (W 3 x H 6 in.)  
**Medium (Dual Direction)** (W 4.5 x H 6.5 in.)  
**Large** (W 9 x H 6.5 in.)
- Allows any Molle/PALS kit to be detachable
- Gives user faster access to Molle/PALS equipped kits
- Mil-Spec hook & loop and heavy duty snaps to ensure kit security
- Designed to be compatible with all of NAR IFAK's, including a NEW horizontal mount option in the medium version (W 4.5 x H 6.5 in.)

**Color Options:**



*Horizontal & Vertical Mounting capability available with the 4.5 x 6.5 in. Panel*

**NAR Quick Detach Panel (3 x 6 in.)**

Color	Item #	NSN#
BLK	80-0893	
COY	80-0894	

**NAR Dual Direction Quick Detach Panel (4.5 x 6.5 in.)**

Color	Item #	NSN#
BLK	80-0895	
COY	80-0896	

**NAR Quick Detach Panel (9 x 6.5 in.)**

Color	Item #	NSN#
BLK	80-0897	
COY	80-0898	



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**

# NAR QUICK DETACH PANELS *3 size options*

REV103117

## NAR Quick Detach Panel (3 in. x 6 in.)



### Features:

- NAR's smallest detach panel measuring (W 3 x H 6 in.)
- Utilizes only 3 rows by 2 columns on user vest/bag, vertical mount configuration
- Designed specifically for NAR's most compact M-FAK kits
- Mil-Spec hook & loop and 3 heavy duty snaps to ensure kit security

### Dimensions:

- W 3 in. x H 6 in. x D 0.5 in.
- Weight: 2.8 oz

## NAR Dual Direction Quick Detach Panel (4.5 in. x 6.5 in.)



### Features:

- NAR's medium size detach panel with dual direction mount measuring (W 4.5 x H 6.5 in.)
- Utilizes 3 rows by 3 columns vertical mount or 2 rows by 4 columns horizontal mount
- Designed for compatibility with all NAR IFAKs such as horizontal mount for M-FAK and vertical or horizontal mount of the TORK
- Mil-Spec hook & loop and 3 heavy duty snaps to ensure kit security

### Dimensions:

- W 4.5 in. x H 6.5 in. x D 0.5 in.
- Weight: 4.5 oz

*Horizontal & Vertical Mounting capability available with the 4.5 x 6.5 in. Panel*



## NAR Quick Detach Panel (9 in. x 6.5 in.)

### Features:

- NAR's large detach panel measuring (W 9 x H 6.5 in.)
- Utilizes 3 rows by 7 columns on user vest/bag, vertical mount configuration
- Designed specifically for use with NAR 4 Chest Pouch & Emergency Evacuation Litter (EEL)
- Mil-Spec hook & loop and 4 heavy duty snaps to ensure kit security

### Dimensions:

- W 9 in. x H 6.5 in. x D 0.5 in.
- Weight: 5.6 oz



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH  **MISSION®**