

QUIKLOT® EMS HEMOSTATIC BANDAGES

REV060823

The QuikClot® EMS Starter Pack provides the tools you need to help promote faster bleeding control on bleeding injuries such as Smaller Surface (Nuisance) Bleedings, Deeper and/or Large Surface Bleedings, to Massive Hemorrhage Penetrating Wound Bleedings encountered in the field, this configuration will address the right product for the right bleeding injury.

Starter Pack Includes:

- QuikClot® EMS 4x4 Dressings (x5)
- QuikClot® EMS Rolled Gauze (x1)
- QuikClot® Combat Gauze® LE (x1)

Features:

- **Controls Bleeding:** Promotes clotting within minutes, allowing responders to control bleeding in an efficient and timely manner.
- **Easy To Use:** QuikClot® EMS Hemostatic Gauze uses the same placement technique as standard gauze
- **Heat Free:** QuikClot® products produce NO Heat

4x4 Dressing Dimensions:

- Packaged - H 5.5 in. x W 5 in. x D 0.2 in.
- Gauze – W 4 in.x L 4 in.
- Weight: 0.8 oz

Combat Gauze LE Dimensions:

- Packaged - H 5 in. x W 5.5 in. x D 0.625 in.
- Gauze – W 3 in.x L 4 yd (12 ft)
- Weight: 0.8 oz



QuikClot® EMS Hemostatic Bandages

Option	Item#	NSN#
EMS Starter Pack	30-0284	
EMS 4x4	30-0095	
Combat Gauze® LE	30-0109	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®