

DRY STERILE BURN DRESSING *Cravat XL*

REV021122

Non-Adherent Dressing...

The Dry Sterile Burn Dressing Cravat XL from North American Rescue is a latex free, non adherent dressing made of woven polyester with two (2) sewn reinforcement seams and high strength polyester lock stitch.

Its generous 45 in. x 45 in. x 63 in. size can be used as a triangular bandage, swathe, sterile field dressing and sling, among other uses.

Vacuum packaged in a 3 mil nylon poly barrier bag, it includes 2 (two) No. 3 closed safety pins and rubber band individually packaged in their own 2 mil nylon poly bag.

Features:

- Dry sterile, non adherent burn dressing and bandage
- Made of woven polyester with two (2) sewn reinforcement seams and high strength polyester lock stitch
- Latex Free, single use
- Vacuum packed in a 3 mil nylon poly barrier bag
- Includes 2 (two) No. 3 closed safety pins and rubber band individually packaged in their own 2 mil nylon poly bag
- Ideal for use as a triangular bandage, swathe, sterile field dressing and sling

Dimensions:

- Bandage: 45 in. x 45 in. x 63 in.
- Vacuum Packaged:
W 4.5 in. x H 4 in. x D 0.5 in.
- Weight: 1 oz

Dry Sterile Burn Dressing Cravat XL

ITEM #	NSN #
30-0162	6510-01-587-6579



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**