

AMBU SPUR II - PEDIATRIC

REV061323

AMBU® SPUR® II provides users with exceptional tactile and visual feedback during resuscitation. The bag is highly responsive, with minimal mechanical resistance. Furthermore the characteristic design provides optimum stroke volume with perfect recoil.

The SPUR® II has a stroke volume of 450 ml and a resuscitator volume of approximately 670 ml.

Features:

- Single-shutter valve system for reliable functionality
- Integrated handle for user comfort and uniform compression
- Thin and responsive bag with minimal mechanical resistance
- SafeGrip™ surface for secure handling in stressful environments
- M-Port allows sidestream measuring of et CO2 or quick medication delivery without disconnecting the AMBU® SPUR® II from the ET tube
- Easy attachment of manometer and PEEP valve
- MR Conditional

Dimensions:

- L 223 mm. x W 99 mm.
- Weight: 145 g.



AMBU® SPUR® II - Pediatric

ITEM #	NSN #
10-0055	



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®