

RAPID RESPONSE KIT *Law Enforcement edition*

REV012820



Rapid Response Kits can be modified upon special request.

Contact a North American Rescue sales representative for more details.

Treat Two Injured Patients Quickly...

Modular, personalized organization of mission essential gear AND medical equipment needed in response to an Active Shooter scenario. Exterior of the rugged nylon shoulder bag features adjustable, padded strap with quick release buckle.

Multiple storage capabilities including loop-tape/MOLLE attachment system, 2 extra pockets & elastic storage loops with color-coded labels to accommodate up to 8 chem lights (not included). Interior compartment features 3 fixed mesh pockets and repositionable elastic loop panel for secure equipment storage. Available in 500D CORDURA® nylon.

The Law Enforcement edition of the Rapid Response Kit is equipped for self aid/buddy aid & configured for treating up to 2 injured patients. MOLLE'd on the exterior are 2 C-A-T® tourniquets with holders, a double mag pouch to accommodate 2 pistol/2 rifle magazines(magazines not included) and roomy nylon pouch to store additional equipment.

Features:

- Modular, personalized organization of mission essential gear AND medical equipment
- Rugged nylon shoulder bag features adjustable, padded strap with quick release buckle
- Multiple storage capabilities including loop-tape/MOLLE attachment system, 2 extra pockets & elastic storage loops with color-coded labels to accommodate up to 8 chem lights (not included)
- Also available as a bag only with shoulder bag, 2 TQ holders, 1 IPOK® Holder, and 1 Double Mag Pouch (no contents included)

Rapid Response Kit - LE Edition

Color	Item #	NSN#
BLK	85-0715	

Rapid Response Bag only - LE Edition

Color	Item #	NSN#
BLK	80-0414	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

RAPID RESPONSE KIT *Law Enforcement edition*

REV012820



QTY
2 C-A-T®
Tourniquet



QTY
2 C-A-T®
TQ Holder



QTY
2 ETD™ 6 in.
Dressing



QTY
2 Gauze,
Compressed



QTY
2 HyFin®
Vent. Twin Pk



PAIR
2 Nitrile
Gloves, Lg



QTY
1 Trauma
Shears, Lg



QTY
1 Double
Mag Pouch



QTY
1 Nylon
Equip. Bag

Kit Contents:

- 2 x Tourniquet Holder
- 1 x Double Mag Pouch
- 1 x IPOK® Holder
- 2 x C-A-T® Tourniquet, Black
- 2 x Emergency Trauma Dressing, 6 in.
- 2 x NAR Compressed Gauze
- 2 x HyFin® Vent Chest Seal, Twin Pack
- 2 x pair Black Talon® Nitrile Glove, Large
- 1 x NAR Trauma Shears, 7.25 in.

Dimensions:

- Kit Closed: H 12 in. x W 12 in. x D 6 in.
- Adjustable shoulder strap:
28 in. to 50 in.
- Weight (with medical supplies):
3 lb 12.5 oz



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®