

# QUANTUM TRAINER CONTROLLER

REV061423

The Quantum Training Controller is specifically designed for use with the Quantum Trainer System ONLY. It provides students with the simulated experience of using the Quantum system for warming IV fluids and blood without actually heating. Lights on the Quantum Trainer Controller unit behave like the real unit to familiarize users with the LED indicators when in real world use. The lights flash blue during the 24 second warm up time and then shift to solid green indicating that fluids would have achieved the proper temperature.

## **Specifications:**

- Allows medical trainers to thoroughly instruct personnel on the proper usage of the Quantum™ Blood and Fluid Warming System. NOT FOR CLINICAL USE. The Quantum™ Trainer System is designed for instructional and training purposes ONLY.

## **Dimensions:**

- Trainer Controller Body: W 1.25 in. x H 3 in. x D 0.75 in.
- Trainer Controller Total Length: 51.25 in.
- Weight (Trainer Controller only): 2.5 oz
- Packaged Weight: 5 oz.

## **Quantum™ Trainer Controller**

ITEM #	NSN #
35-0027	----



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**