

QUANTUM TRAINER TUBING SET

REV 061523

The Quantum Trainer Tubing Set is a component for use with the Quantum Trainer System ONLY. It is designed for non-warming training purposes and does not contain any heating elements. Designed to provide users with the necessary equipment to confidently practice the setup process and connection techniques, simulating a real-life scenario. The Quantum Trainer System allows for users to gain hands-on experience and develop proficiency in utilizing the actual Quantum™ Blood and Fluid Warming System. By mastering the setup and operation procedures through practical training, healthcare professionals can ensure efficient and effective fluid warming practices when using the actual Quantum™ Blood and Fluid Warming System in clinical settings.

Specifications:

- Allows medical trainers to thoroughly instruct personnel on the proper usage of the Quantum™ Blood and Fluid Warming System. NOT FOR CLINICAL USE. The Quantum™ Trainer System is designed for instructional and training purposes ONLY.

Kit Contents:


- 1 x Y-Type Dual Spike
- 1 x 200 micron filter
- 1 x Male Luer Lock adapter
- 2 x Luer activated injection sites

Dimensions:

- Trainer Tubing package (each): 6.5 in. x 11 in.
- Tubing Length: 85 in.
- Trainer Tubing Weight (each): 2.7 oz

Quantum™ Tubing Trainer

| ITEM # | NSN # |
|---------|-------|
| 35-0028 | ---- |



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**