

TCCC/TECC Modular Training System | Module 3: Extraction

REV111220



QTY 1 Rat™ Strap
 QTY 1 Dragon Handle System
 QTY 1 Hasty Harness
 QTY 1 Non-Locking Carabiner



QTY 1 Extraction Device, TED
 QTY 1 Litter, EEL
 QTY 1 Talon 90C Litter

Module Three – Extraction: Provides various high-threat extraction devices and casualty evacuation systems to train students in different options to move casualties during the all three TCCC Phases of Care. Packaged in a self-contained, wheeled, rigid case that allows for easy transport and storage, this kit includes High Threat Extraction drag devices, versatile casualty movement Hasty Harness, poleless/non-rigid litter option and the collapsible quad-fold rigid evacuation platform standard – Talon II 90C Litter.

Hard Case Specifications:

- Retractable extension handle
- 4 strong polyurethane wheels with stainless steel bearings
- Easy-open double throw latches
- Open cell core with solid wall design – strong & lightweight
- O-ring seal
- Automatic pressure equalization valve
- Fold down handles
- Stainless steel hardware and padlock protectors
- Desert tan color

Dimensions:

- Hard Storage Case:
 External – L 30.75 in. x W 20.5 in. x D 11.6 in.
 Internal – L 28.5 in. x W 17.5 in. x D 10.5 in.
- Weight: Module Three - 45 lb



Module Three Contents:

- 1 x Hard Storage Case
- 1 x Talon II® 90C Collapsible Handle Litter
- 1 x EEL, Emergency Evacuation Litter- BLK
- 1 x TED, Tactical Extraction Device- BLK
- 1 x Hasty Harness - Blk
- 1 x Dragon Handle System
- 1 x RAT Strap, Rescue Assault Tether
- 1 x Shield Carabiner

TCCC Modular Training System Module Three: Extraction

Item #	NSN#
85-0521	



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A
MISSION®