

# LATERAL CANTHOTOMY SET

REV061623

The ability to perform a lateral canthotomy and cantholysis is an essential skill on the battlefield in that it can be a vision-saving technique in the setting of decreased perfusion to the optic nerve or globe as a result of increased orbital pressure. This Lateral Canthotomy Set contains everything a provider would need to perform these procedures.

## Specifications:

- Compact kit allowing for ease of use in the field
- Tools to decrease orbital pressure

## Kit Contents:

- 1 x Eye Shield, Polycarbonate - PES - Orange
- 4 x Gauze, Pad - 4 in. x 4 in. (Pk of 2)
- 1 x Halstead Mosquito Forceps - Straight 5 in.
- 1 x Iris Scissors, Straight - 4.5 in.
- 1 x Metzenbaum Scissors, Straight 5.5 in.
- 1 x Syringe, Tuberculin with Needle 1ml (27 ga x 1/2 in.)
- 1 x Adson Tissue Forceps, 1 x 2 - Sterile

## Dimensions:

- L 8 in. x W 4.5 in. x D 2 in.
- Weight: 5.5 oz

## Lateral Canthotomy Set



ITEM #	NSN #
83-0012	--



## Kit Contents:



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
 WITH  **MISSION®**