

EZ-IO® TRAINING NEEDLES *Red Training Versions*

REV100620

Red Colored Training Needle Hub Versions...

EZ-IO® Power Driver Training Needles provide medical professionals with the the ability for complete training on the EZ-IO® Power Driver Trainer Model.

EZ-IO® provides the medical professional immediate vascular access to the central circulation within seconds, delivering medications, intravenous fluids and blood products to adult and pediatric patients. With a specially designed cutting IO needle, manual driver, or small power driver, the EZ-IO® allows the clinician complete control — avoiding the use of force. EZ-IO® provides rapid, smooth entry into the bone's medullary cavity, creating an immediate conduit to the central circulation.

The EZ-IO® Intraosseous Infusion System is the leading choice for intraosseous access across multiple health care settings, and is used by 90 percent of US advanced life support ambulances, over half of US Emergency Departments, and the US Military.

Unlike older IO technology, the EZ-IO®'s patented needle tip has four asymmetrical bevels on the stylet, and two cutting tips on the catheter, providing a smooth and accurate insertion. Black indicator lines are located at 5 mm increments from the hub for assessing insertion depth. With EZ-IO®, insertion depth is completely controlled by the clinician.

Three needle sizes accommodate intraosseous access in patients 3 kg and larger.

The red coloring of the Trainer Version Needles provides an immediate visual clue that they are only intended for training purposes. The trainer needle versions feature the same weight and dimensions as the tactical use versions so trainees and medical personnel can familiarize themselves with the proper use and techniques required for in-field applications.

Available in three sizes:

15 mm Needle Set (Item #30-0081) Designed for use in patients weighing 3-39 kilograms, and for patients with minimal tissue over insertion sites.

25 mm Needle Set (Item #30-0083) Designed for any patient weighing more than 39 kg or for patients who have too much tissue over the insertion site for the 15 mm Needle Set to be used.

45 mm Needle Set (Item #30-0085) Designed for patients weighing more than 39 kg who have excessive tissue over the targeted insertion site (e.g. edema, large musculature, obesity). The 45 mm Needle Set is ideal for the humerus site in patients over 39kg

EZ-IO Products are restricted to Federal Purchases only. Please contact Customer Service for details.

**Caution: Federal (USA) law restricts this device to sale by or on the order of a practitioner licensed to use or order use of the device or Federal (USA) law prohibits dispensing without a prescription.*



Dimensions:

- Packaged: H 4.75 in. x W 4.25 in x 0.7 in.
- Weight: 1 oz

EZ-IO® Training Needles - Red Hub

Size	Item #	NSN#
15 mm Pediatric	30-0082	
25 mm Adult	30-0084	
45 mm Adult Large	30-0086	

CE



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®