

# EZ-IO® G3 POWER DRIVER *Trainer*

REV113018

*EZ-IO Products are restricted to Federal Purchases only.  
Please contact Customer Service for details.*

*\*Caution: Federal (USA) law restricts this device to sale by or on the order of a practitioner licensed to use or order use of the device or Federal (USA) law prohibits dispensing without a prescription.*



Molded Hard Side Carrying Case sold separately

CE

## **EZ-IO® G3 Power Driver Trainer**

Unit	Item #	NSN#
Driver	30-0076	
Case	30-0077	

### **Specifications:**

#### **Power Driver:**

- H 4.5 in. x W 6.5 in. x D 2.5 in.
- Weight: 15 oz

#### **Hard Case:**

- L 8.5 in. x W 6 in. x D 2.25 in.
- Weight: 0.65 lb.

### **Red Color Trainer Version...**

EZ-IO® Power Driver Trainer provides medical professionals with the training needed to be able to effectively and efficiently use the EZ-IO® G3 Power Driver in the field.

The EZ-IO® G3 Power Driver provides immediate vascular access to the central circulation within seconds, delivering medications, intravenous fluids and blood products to adult and pediatric patients. With a specially designed cutting IO needle, manual driver, or small power driver, the EZ-IO® allows the clinician complete control — avoiding the use of force. EZ-IO® provides rapid, smooth entry into the bone's medullary cavity, creating an immediate conduit to the central circulation.

The EZ-IO® G3 Power Driver – applicable to all levels of care – has 15mm, 25mm and 45mm IO needle options. It can provide vascular access to peripheral intraosseous access sites (humerus and tibia) in less than 10 seconds for both adult and pediatric patients. The weight and compact design of the G3 IO driver allows medical personnel to carry enough components for multiple insertions with a fraction of the weight of a comparable number of other IO products.

The red coloring of the Trainer Version provides an immediate visual clue that it is only intended for training purposes. The trainer version features the same weight and dimensions so trainees and medical personnel can familiarize themselves with the proper use and techniques required for in-field applications.

### **Features:**

- Battery-powered revolutionary technology
- Provides access within 10 seconds
- Multiple insertion sites with three needle sets
- Standard Luer-Lock catheter for fast, easy removal
- Trigger guard prevents accidental activation of the device when stored in kits
- Lightweight, compact design in a distinct Red Training color
- Molded Hard Side Carrying Case (sold separately - Item # 30-0077)



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**