


FOX EYE SHIELD

REV062123



Convex Rigid Eye Shield...

The Rigid Eye Shield meets current military and civilian pre-hospital guidelines for the treatment of traumatic eye injury in the field.

This convex rigid eye shield allows the user to protect the injured eye and blocks external pressure from being applied directly to the injury, preventing further damage.

It is made of aluminum for rigidity, yet contains many apertures to allow for air circulation and to reduce weight.

Features:

- Combination Eye Garter Shield Cover and Fox Shield; includes E5692 Pack of 12 eye shields
- Lightweight aluminum
- Perforated to allow air circulation
- Accommodates either eye
- Meets CoTCCC treatment protocols for penetrating eye trauma

Dimensions:

- Dimensions: 2.938 in. x 2.312 in.
- Individual Weight: 0.2 oz (5 grams)

Fox Eye Shield

CE

OPTION	QTY	ITEM #	NSN #
No Garter	Box of 12	30-0144	6515-01-449-1016



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®