

PVC BLEEDING CONTROL PATCH

REV062723

The PVC Bleeding Control Patch from NAR is the perfect accessory for labeling your lifesaving medical gear. Originally designed for use on the Public Access Bleeding Control Kit, this durable PVC rubber patch features 3D easy-to-read letters and hook field backing that can be used on any surface with an area for hook-&-loop attachment. An additional loop field lined sheet is included to protect the hook back when patch is not in use.

Features:

- Soft flexible PVC rubber patch for labeling your lifesaving medical gear
- Ideal for the Public Access Bleeding Control Kit or any surface with an area for hook-&-loop attachment
- Backed with hook material for ease of attachment and removal
- Easy-to-read 3D letters in a convenient 4.375 in. x 1.75 in. size
- Additional loop field included to protect hook back when patch is not in use

Dimensions:

- L 4.375 in. x H 1.75 in. x D 0.125 in.
- Weight: 0.275 oz



PVC Bleeding Control Patch

ITEM #

ZZ-0796



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®