

# RIG SERIES Assault Litter Carrier

REV101818



NOTE: Bag only.  
Litter, HRS, Heating  
Element & Casualty  
Equipment Bags all  
sold separately.



## New & Enhanced Assault Litter Carrier...

New enhanced version of our standard Assault Litter Carrier with additional securing strap over the top of the litter to hold it in place even during drastic movement and FRIES (Fast Rope Insertion Extraction System) operations.

Upgraded shoulder straps provide comfort for long-term wear. Includes pocket in the top flap for storage of Heat Reflective Shell (HRS) & pocket in the back panel for the Heating Element & Casualty Equipment Bag (all items sold separately).

Three additional handles facilitate packing & unloading of the litter. Equipment attaching webbing (MOLLE) along the outside of the carrier allows attachment of additional mission dependent equipment.

Available in 500D CORDURA® Solution-Dyed Nylon and GhillieTEX™ hardware for IR signature reduction.

Accommodates all NAR quad-fold litter platforms.

### Features:

- Litter Carrier with additional securing strap over the top of the litter to hold it in place even during drastic movement and FRIES (Fast Rope Insertion Extraction System) operations
- Upgraded shoulder straps provide comfort for long-term wear
- Includes pocket in the top flap for storage of Heat Reflective Shell (HRS) & pocket in the back panel for the Heating Element & Casualty Equipment Bag (all items sold separately).
- Three additional handles facilitate packing & unloading of the litter

### Specifications:

- Dimensions:  
H 22 in. x W 9 in. x D 9 in.
- Weight: 3 lb

### RIG Series Carrier ONLY

\*NATO NSN

Color	Item#	NSN#
COY	60-0056	6530-01-643-9138*
MTC	60-0057	6530-01-647-8571*



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**