

T.A.L.O.N.™ MANUAL TRAINING NEEDLE SET (BOX OF 6)

BY EZ-IO®

REV062723

Manual I/O Insertion (Training)...

Non-sterile device for the training of Manual I/O insertion at all 7 extremity sites: Sternum, Bilateral proximal humerus, Bilateral proximal tibia, and Bilateral distal tibia. Includes 15 gauge, 304 stainless steel Needle/Catheter, Sharps Protector, EZ-Connect® Extension Set, I/O Wristband, Quick Reference Card & Sternal Locator. Also comes with re-set tool for removal and re-use.

NOTE: Non-Sterile, for training purposes only.

*** EZ-IO® Products are restricted to Federal Purchases only. ***

*** Please contact Customer Service for details. ***

Features:

- Sold in Box of 6
- Allows for training in manual I/O insertion at all 7 extremity sites: sternum, bilateral proximal humerus, bilateral proximal tibia, and bilateral distal tibia
- Includes 15 gauge, 304 stainless steel Needle/Catheter, Sharps Protector, EZ-Connect® Extension Set, I/O Wristband, Quick Reference Card & Sternal Locator
- Also comes with re-set tool for removal and re-use
- NOTE: Non-Sterile, for training purposes only

Specifications:

- Packaged: L 10 in. x W 8.5 in.
- Weight: 6.5 oz



T.A.L.O.N.™ Manual Training Needle Set

| ITEM # | NSN # |
|---------|-------|
| 30-0133 | ---- |



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**