

IV POLE *Adjustable*

REV062722



Features:

- Essential, lightweight litter accessory designed to fit the pre-molded holes of our Talon® and Raven® framed litters
- Easy to transport and fully adjustable, from a minimum clearance of 18.5 in. expanding to a height of 32 in. when deployed

Efficient and Reliable Soution...

North American Rescue's Adjustable IV Pole is exclusively designed to expand the treatment capabilities of our family of evacuation platforms by fitting the pre-molded holes of all NAR framed litter platforms including the Talon®, Raven®, T3, Medevac & OSL litters.

This essential, lightweight accessory is easy to transport and fully adjustable, from a minimum clearance of 18.5 in. expanding to a height of 32 in. when deployed.

Exclusively designed to fit the pre-molded holes of our framed litters and weighs less than 5 ounces.

Specifications:

- Retracted: H 18.5 in. x W 8.86 in. x D 0.5 in.
- Fully Extended Height: 32 in.
- Weight: 4.9 oz.

NOTE: All accessories are sold separately. Litters are not included.

Adjustable IV Pole

Item#	NSN#
60-0008	6530-01-504-9054
CE 60-0008-EU	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®