

NAR HEAD IMMOBILIZER

REV120518



*Easy attachment to
spine board
(sold separately)*

Easy to use, compact and extremely effective...

Constructed of soft, closed-cell foam with hook and loop padded straps, the NAR Head Immobilizer provides quick and effective cervical restriction.

Ears and head are accessible without removing the device for ongoing patient assessment.

Features:

- Provides quick and effective cervical restriction
- Soft, closed-cell foam with hook and loop padded straps
- Ears and head are accessible without removing the device

Specifications:

- H 10 in. x W 16 in. x D 5.75 in.
- Weight: 2 lb 2.8 oz
- SEL#: 09ME-04-SPIN

NAR Head Immobilizer  

Item#

NSN#

50-0020

6530-01-619-1777



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**