

I&A CASE *Injury & Ailment Case*

REV062923

Lid's Clear Pockets And Hook-And-Loop Fasteners Allow Secure Storage and Quick Inventory...

Engineering excellence is a hallmark of North American Rescue. With our Injury and Ailment Case (I&A), engineering excellence is evident in every detail. Designed for healthcare professionals on the move, this pelican case is a tribute to the experts who developed it as well as to the experts who use it.

Engineering excellence prompted a lid organizer that opens double for double capacity. The lid's clear pockets with hook-and-loop attachments allow quick inventory and efficient storage of medications. The lids have other advantages, too. A polymer O-ring creates a seal that is dust-proof and waterproof. The rubberized sure-grip handle is reliable in any weather. Engineering excellence made this med case exceptionally secure. The case is lockable, of course, but that's just the beginning.

Stainless steel reinforced padlock protectors provide a second-layer of defense against cutting and theft. The med box also features double-throw latches which are engineered to tolerate up to 400 pounds of pressure, yet open with a light pull. The case's standard automatic Gore-TEX pressure equalization valve not only stops moisture from entering the case, but also prevents vacuum-lock. These medicine cases open easily at any altitude, but only when they are supposed to open.

The main body of the case is equipped with a fully customizable insert covered with 210 Denier Oxford nylon. The insert's padded dividers (made of half-inch PE foam) have adjustable walls for customizing compartments to any size.

Flexible hook-and-loop attachment points accommodate the most demanding requirements. Adjustable walls plus convenient attachment equals extra protection and superior inventory management. Our medicine cases are large enough to handle the job, and small enough to transport with ease.



NOTE: Medical equipment is not included.

Injury & Ailment Case

ITEM #	NSN #
ZZ-0222	



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**

Features:

- SEL#: 09ME-01-BAGM
- Lid's clear pockets and hook-and-loop fasteners allow secure storage and quick inventory
- Insert cover is 210 Denier Oxford Nylon
- Insert's padding and flexible fasteners support customization of the case's contents
- Lid's polymer O-ring gives a dust- and water-proof seal
- All-weather, sure-grip handle
- Reinforced padlock protectors

Dimensions:

- Exterior: H 8.5 in. x W 20.5 in. x D 17 in.
- Interior: H 7.6 in. x W 18.4 in. x D 14 in..
- Weight: 13 lb 12 oz

Superbly engineered medical equipment organization system for healthcare professionals on the move.



NOTE: Medical equipment is not included.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A 
MISSION®