

# COURTESY MASK

REV070523

The NAR Courtesy Mask is a comfortable cloth mask with elastic bands that loop over the ears. The fabric is pleated so that it easily conforms to the face, wrapping over the nose and under the chin for maximum coverage.

The 2-ply fabric is 100% cotton muslin that is washable and re-usable. The knitted elastic bands are sewn into loops that go over the ears, applying sufficient tension to secure the fabric firmly, while offering a comfortable fit to the face.

Care: Hand wash (recommended) or Machine Wash Cold. Hang dry. Do not heat, as this will damage the elastic straps, causing them to lose their elasticity. Use non-chlorine bleach to avoid damaging the color.

## Features:

- Comfortable 100% cotton muslin
- 2-Ply pleated face cover design
- Knitted elastic bands loop over ears for secure fit
- One size fits most

## Specifications:

- Mask: 7.25 in. x 3.75 in.
- Bands: 8.00 in. x 0.25 in.
- Weight: 0.01 lb



## Courtesy Mask

ITEM#	UNIT
ZZ-1316	Single



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**