

OPHTHALMOSCOPE *Heine Mini 3000*

MED/SURGE



Designed for examining the interior structures of the eye, especially the retina, in a mini, pocket-size format.

High-Quality Optic Tool...

The Heine® 3000 Ophthalmoscope is designed for examining the interior structures of the eye, especially the retina, in a mini, pocket-size format.

This critical high quality optic tool offers all of the important features of a high-end ophthalmoscope combined in a compact, durable, portable configuration.

Shockproof, sturdy, non-slip and easy to maintain.

Features:

- Enhanced XHL Xenon Halogen Technology; concentrated white light for perfect illumination & crisp retinal detail
- -20 to +20 with 18 lenses; -/+1 to -4/+4 in single steps
- 5 apertures for full diagnostic capability
- 2-piece head & handle with chrome finish upper section/ refined plastic; Shockproof, sturdy, non-slip
- Easy to maintain, flexible
- Integral pocket clip with auto off
- Takes 2 AA batteries (included)

Dimensions:

- H 6.75 in. x W 1.25 in. x D 0.9 in.
- Weight: 3.7 oz

Heine Mini 3000 Ophthalmoscope

Item#	NSN#
20-0007	6540-01-491-6651

High quality optics for examining the interior structures of the eye in a mini, pocket size.

PRODUCTS WITH A
MISSION®



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277