

# EYESHIELDS WITH SELF ADHESIVE *Pack of 100*

REV070623

Personal protective equipment (PPE) is crucial in order to prevent the spread of infectious viruses, bacterial infections and other diseases. This lightweight transparent polyester eye shield includes a peelable liner with 2 self adhesive strips designed to adhere to any face mask or respirator (sold separately). This single use eye shield is manufactured in a clean ISO approved facility and is designed to protect the face from splashing or spraying of infectious or hazardous substances. The tested unique design means for a comfortable fit when worn for prolonged periods of time. \*Eye Shield for use in conjunction with mask or respirator (sold separately) and sold in a pack of 100.

## Features:

- 3M™ 6.9Mil Anti-Static
- Transparent (optically clear)
- Non-Hemolytic
- Low Autofluorescence
- Non-toxic to mammalian cells
- No significant outgassing <400°C
- Single use item
- Material: Polyester
- Made in the USA

## Application Technique:

- Peel liner from the 2 adhesive straps
- Apply eye shield to mask with adhesive facing away from you applied to the outside of the mask
- Apply firm pressure for ultimate bond

## Specifications:

- Dimensions: W 10 in. x H 4.75 in.
- Dimensions (packaged): W 12 in. x H 7 in.
- Weight: 0.32 oz (each) / 2 lb (pack of 100)



## *Eyeshields with Self Adhesive*

ITEM#	UNIT
ZZ-1351	Pack of 100



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH  
**MISSION®**

