

DRAGON HANDLE SYSTEM

REV041522

Handle Assisted High Threat Extraction...

This innovative extraction system is constructed of 1 in. Tech Tape with a tensile strength of 4000 lb, outperforming 1 in. military spec webbing not only in strength but also in flexibility, knotability and durability.

Two weather-proof flexible handles offer additional grip and comfort to the rescuer. The floating black bull ring provides adjustable attachment point configurations.

Also includes a Black Omega, non-locking "D" carabiner rated for 31kN (6968.8 lb)

Features:

- Designed for linear personnel recovery drags
- Decreases time on the "X"
- Exponentially faster drag time with less physical exertion
- Keeps victim low profile, decreasing the likelihood of additional injuries
- Head first or Feet first drag options
- Patent: 8,016,335

Dimensions:

- Packaged: H 6.3 in. x W 4.1 in. x D 2.25 in.
- Extended Length: L 42 in.
- Minimum Break Strength: 1443 lb
- Weight: 11.3 oz

Warning: Not designed for hoisting or lowering operations

All items must be visually inspected for any deterioration of stitching and webbing both before and after each use.



Dragon Handle System



Item#	NSN#
90-0004	6530-01-598-9595



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®