

LACERATION TRAYS

Facial Laceration, Minor Laceration & Toenail Removal

REV010323

Cost effective, sterile, single use kits with a 3 year shelf life...


Save time and eliminate waste by having just the right components at your finger tips to provide patients with the treatment they need faster.

Sterile, ready-to-use disposable kits provide the essential latex-free supplies and instruments needed for minor surgical procedures commonly performed in the field hospital.

Offering extra protection against cross infection, these compact convenient kits help increase efficiency while reducing set up time and cost.

Facial Laceration Tray

ITEM #	NSN#
ZZ-0683	6515-01-661-1802



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

Features:

- Sterile
- Easy Open Pouch
- Cost Effective

STERILE



Specifications:

- Packaged: L 9 in. x W 5 in. x D 2.75 in.
- Weight: 12 oz

**Customizable Procedure Trays/Kits are also available upon request. To create a kit to meet your exact requirements and specific needs, please contact North American Rescue for additional information.*



Facial Tray Kit Contents:

- 1 x Kinguard One Step (30 in. x 30 in.)
- 1 x Medicine Cup (60cc)
- 1 x Scalpel, Safety #15
- 1 x Halstead Mosquito Forceps, ST Satin
- 1 x Halstead Mosquito Forceps, CRV Satin
- 1 x Webster Needle Holder, Satin (5in.)
- 1 x Adson Forceps, FP, Satin
- 1 x Iris Scissors, ST, Satin
- 1 x Iris Scissors, CRV, Satin (4.5 in.)
- 1 x Frazier Skin Hook, Sharp (5 in.)
- 10 x Gauze Pads 4 x 4 in. 12 ply
- 2 x Gauze Pads 2 x 2 in. 12 ply
- 1 x Fenestrated Drape with tape (18 in. x 26 in.)
- 2 x Towel, OR Blue (17 in. x 25 in.)

PRODUCTS

WITH A

MISSION®

LACERATION TRAYS

Facial Laceration, Minor Laceration & Toenail Removal

REV010323

Toenail Removal Kit

ITEM #

ZZ-0685



Toenail Removal Kit Contents:

- 1 x Kimguard One Step (30 in. x 30 in.)
- 1 x Medicine Cup (60cc)
- 1 x Probe, with Eye (5 in.)
- 1 x Kelly Hemostat, ST (5.5 in.)
- 1 x Adson Forceps, 1x2 (4.75 in.)
- 1 x O.R. Scissors, SB ST (5.5 in.)
- 12 x Gauze Pads 4 x 4 In. 12 ply
- 1 x Fenestrated Drape with tape (18 in. x 26 in.)
- 2x Towel, Absorb. White (17 in. x 20 in.)




Minor Laceration Kit

ITEM #

NSN#

ZZ-0684

6515-01-661-1786



Minor Laceration Kit Contents:

- 1 x Scalpel, Safety #15
- 1 x Halstead Curved Mosquito Forcep (5 in.)
- 1 x Adson Forceps FP, 1x2 (4.75 in.)
- 1 x Iris Scissors, CRV (4.5 in.)
- 1 x Webster Needle Holder
- 2 x Gauze Pads 2 x 2 in. 12 ply
- 10 x Gauze Pads 4 x 4 in. 12 ply
- 1 x Wrap (30 in. x 30 in.)
- 1 x Fenestrated Drape with tape (18 in. x 26 in.)
- 2 x Towel, Absorb. White 17" x 20"
- 1 x Medicine Cup (2 oz)



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®