

# TOMMANIKIN® RIGHT ARM WITH SHRAPNEL INJURIES

REV071023

## Active Bleeding Simulator

The TOMManikin® Right Arm with Shrapnel Injury features widespread penetrating shrapnel injuries, severe soft tissue damage, and bleeding from several locations. This training simulator is tourniquet responsive, this allows the users to train and utilize the Active Bleed feature of the training simulator. Constructed with 100% silicone and a durable metal bone structure, the Right Arm with Shrapnel Injuries is interchangeable with all TOMManikin® and TAMiKin® simulators.

## Specifications:

- Tourniquet Responsive
- Exposed tissue and muscle layers
- Integrated bleeding system
- Compatible with all TOMManikin® and TAMiKin® simulators

## Dimensions:

- L 24 in. x W 12 in. x H 6 in.
- Weight: 10 lbs.

## TOMManikin® Rt. Arm w. Shrapnel

ITEM #	NSN #
93-0085	----



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A **MISSION**

