

TOMMANIKIN® LEFT ARM UNINJURED

REV071023

An Armload of Features

The TOMManikin® uninjured left arm features radial pulse, antecubital IV access (with flash), and humeral IO access (procedural only). The simulator's pulse rate, pressure, and regularity controlled via NARS tablet.

Features 100% silicone construction with metal bone internals. Interchangeable with all TOMManikin® and TAMiKin® simulators.

Specifications:


- Features radial pulse
- Antecubital IV access
- Humeral IO access (procedural only)
- 100% silicone construction

Dimensions:

- L 24 in. x W 12 in. x H 6 in.
- Weight: 10 lbs.

TOMManikin® Left. Arm - Uninjured

ITEM #	NSN #
93-0088	----



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**