

TOMMANIKIN® RIGHT LEG WITH PARTIAL AMPUTATION

REV071023

An Armload of Features

TOMManikin® right leg partial amputation (tourniquet responsive). Features massive bleeding, exposed bone and extensive soft tissue damage.

This partial amputation trainer is interchangeable with all TOMManikin® and TAMlkin® simulators. 100% silicone construction with metal bone structure.



Specifications:

- 100% silicone construction with metal bone
- Interchangeable with all TOMManikin® / TAMlkin®
- Tourniquet Responsive

Dimensions:

- Packaged - L 36 in. x W 18 in. x H 18 in.
- Weight: 26 lbs

TOMManikin® Right Leg - Part. Amp.

ITEM #	NSN #
93-0096	----



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**