

CREDO PROMED™ SERIES FOUR CONTAINER

REV071223



Dimensions:

2 Liter

- L 8 in. x W 9 in. x D 10 in.
- Weight: 7.64 lbs

4 Liter

- L 10 in. x W 10.5 in. x D 10.5 in.
- Weight: 11.45 lbs

8 Liter

- L 11 in. x W 15.5 in. x D 11 in.
- Weight: 7.64 lbs

Proven Payload Protection...

Designed specifically for the dynamic needs of first responders, the Credo ProMed™ product line currently consists of three different sized temperature-controlled portable medical transport bags that thermally protect the integrity of valuable medical supply payloads for 48 - 96 hours.

The outer bag is constructed of highly durable ballistic nylon fabric and the patented Thermal Insulation Chamber (TIC™) panels with phase change material and VIP components contained within, are qualified to consistently protect medical materials such as blood and platelets within 2-8°C.

These hand portable, highly durable and convenient to use soft-sided bags offer: Easy conditioning of the TIC™ system panels which makes the ProMed a simple, nimble, and consistently reliable option. Once the TIC™ system panels are staged and properly placed inside the Vacuum Insulated Panel (VIP) assembly, simply insert payload content and be on your way.

Rugged and sturdy supportive interior side panels, double-stitched seams, reinforced protective corners, durable buckle straps, comfortable padded shoulder strap and a water and soil resistant protective bottom panel. A sophisticated design with a clear view pocket for business essentials and a sturdy elastic back strap designed to slip effortlessly over telescoping handles enabling secure transport.

These containers will compliment the Quantum and LifeFlow to provide a complete blood delivery package.

Features:

- Duration: Up to 72 hours
- Temperature: 2° to 8° C
- Volume: 2L - 8L
- Outer Container: Ballistic nylon fabric with durable web strapping and secure buckles
- Standard Color: Tan

Crêdo ProMed™ Series Four

CAPACITY	ITEM #
2L	35-0011
4L	35-0012
8L	35-0013



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**