

GOLDEN HOUR™ ONE-V CONTAINER PELICAN™ (0.5L)

REV071323

...from anywhere in the world — to where they need it most...

The Golden Hour™ ONE-V Container is an ultra-lite blood transport capable of carrying one half liter of whole blood and keeping the stored blood cool for up to 48 hours, regardless of weather conditions. This container easily transports a single unit of whole blood to any location. This extremely light-weight container will securely and reliably transport blood and platelets to wherever it is needed most.

Utilizing Thermal Isolation Chamber System (TIC™) technology and vacuum insulation panels, the Golden Hour™ ONE-V Container will keep the payload isolated from exposure to extreme temperatures with a tight, protective thermal seal for a pre-determined temperature range and duration. Safely deliver whole blood and blood supplies from anywhere in the world — to where they need it most. The containers are easy to assemble, reusable parcels that are constructed of reusable non-toxic materials.

Features:

- Ultra-light storage and transport of 0.5L of whole blood
- Keep blood cool for up to 48 hours, regardless of weather conditions
- Dismount for up to 12 hours without external power
- Equip to make blood in far forward areas accessible
- Reusable, rugged design for easy cleaning and replacement

Dimensions:

- W 6 in. x H 8 in. x D 5 in.
- Weight: 2 lbs

Pelican™ Golden Hour ONE-V

ITEM #	NSN #
35-0016	----



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**