

GECKO GRIP MULTIPURPOSE FLAT TAPE

REV071323

The NAR® Gecko Grip Multi-Purpose Flat Tape was designed to provide users with a flat pack configuration for 90 inches of 2-inch-wide high strength T-REX Gecko Grip tape. This is ideal for Outdoor and Adventure IFAKs, and vacuum sealed kits that do not need a full roll of tape. Each kit contains 6 strips, each 2" wide x 15" long with easy peel off backing. The tape kit is lightly vacuum sealed to reduce size but also maintain flexibility for various packing configurations and locations in kits. Contained in a clear resealable poly bag.

Specifications:

- 90 in. of 2 in. wide high strength T-REX Gecko Grip Tape
- Ideal for packing in small kits
- Each pack includes 6 tape strips

Dimensions:

- Packaged: H 7.5 in. x W 3 in. x D 0.34 in. (tape tri-folded inside kit)
- Min Pk Size: H 5.25 in. x W 2 in. x D 0.34 in. (tape tri-folded inside kit)
- Strip Size: 2 in. x 15 in. per strip
- Weight: 0.13 lbs (2.13 oz)

Gecko Grip Multipurpose Flat Tape (6 Pack)

ITEM #

ZZ-1523



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**