

LITTER BUCKLE REPLACEMENT

REV072023

State of the art dual adjust model ITW Nexus TSR 200-BLK features patented Slick-Click latching mechanism and contoured profiles for maximum webbing grip while allowing smooth and easy adjustment. One package includes one male and one female component.

This buckle set is compatible with all of NAR's litters (Talon, Raven, T3, MedEvac and OSL) as well as the NAR Rescue Sleeve II. Buckle set is only available in color Black.

Features:

- State of the art dual adjust model ITW Nexus TSR 200-BLK
- Slick-Click latching mechanism
- Contoured profiles for maximum webbing grip
- Smooth and easy adjustment
- Packaged with one male and one female component
- Compatible with all of NAR's litters (Talon, Raven, T3, MedEvac and OSL) as well as the NAR Rescue Sleeve II
- Color: Black

Dimensions:

- Packaged: L 3.75 in. x W 2.375 in. x D 0.5 in.
- Weight: 1 oz



Litter Buckle Replacement

ITEM#	NSN#
ZZ-0032	6530-01-435-7337



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®