

# OSL™ BACKREST

REV070822



## Efficient and Reliable Solution for Oversized Patients...

For appropriate treatment and improved comfort, supine is not always the best position for patients. When casualties can be better served by being placed at a 45° or 90° angle, the OSL Back Rest offers an efficient and reliable solution.

The OSL back rest attaches to the 2 in. poles of the OSL litter system.

The backrest cover is made of a Nomex/Kevlar blend: Nomex for flame-resistance and heat-energy absorption; Kevlar for abrasion-resistance and long life. The fabric can be repeatedly washed or dry-cleaned without losing these critical properties.

*When the patient's condition requires a 45° or 90° angle rely on the OSL Backrest.*

### Specifications:

- H 31.5 in. x W 30.25 in. x D 1.6 in.
- Weight: 2 lb 15 oz

### Features:

- Covered with durable Nomex/Kevlar-blend fabric
- Positions patients at 45° and 90° angles
- Attaches to the 2 in. poles of the OSL™ litter system
- Increases patient comfort
- SEL#: 09ME-01-COTS

*NOTE: All accessories  
are sold separately.  
Litters are not included.*

### OSL Backrest



Item#	NSN#
60-0024	

*Constructed of a lightweight aluminum  
frame with durable Nomex/Kevlar  
blend fabric cover.*



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**

WITH A

**MISSION®**