

# NAR® MINI RESPONDER ETD® EMERGENCY TRAUMA DRESSING

REV072823



***Our most compact compression bandage designed to provide an effective pressure bandage when reduced packaging size is required...***

The NAR® Mini Responder ETD® from North American Rescue® is a flat folded 42 in. elastic bandage with 4 x 6 inch non-adherent pad designed to be easily stored on the person in pant pockets, shirt pocket, trauma plate pouches, etc. It's small, compact packaging is less than an inch thick with easy open Red-Tip Technology® tear notches. The Mini Responder ETD® has both hook & loop as well as a sturdy C-Clip closure to hold the dressing in place. The bandage is sterile unless opened or damaged.

### Features:

- Small, compact packaging
- Flat folded
- Both hook and loop as well as C-Clip securing device
- 4 x 6 inch non-adherent sterile pad
- 42 inch elastic bandage (unstretched)
- Sterile unless package is opened or damaged

### Dimensions:

- Packaged: L 5 in. x W 3.5 in. x D 0.75 in.
- Bandage length (unstretched): 42 in.
- Pad: 4 in. x 6 in.
- Weight: 2.4 oz

### NAR® Mini Responder ETD®

WIDTH	ITEM #	NSN#
4 in.	30-0179	----



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**