

EMERGENCY TRAUMA DRESSING - ETD®

ABDOMINAL/STUMP ETD®

REV072823



Rapid, controlled, effective bandaging for larger wound patterns

The Abdominal/Stump ETD® is an ideal for dressing abdominal eviscerations, traumatic amputations, burns and large pattern wounds. Engineered using data from over sixty-three after action reviews (AAR) of real-world battle dressing applications, the Abdominal/Stump ETD® consists of a resilient elastic wrap equipped with a 12 in. x 12 in. sterile non-adherent pad, durable securing device and innovative Quick-Grip Roll Control configuration. This combination facilitates a rapid, controlled, and extremely effective application process in the most demanding of conditions.

The ETD® can also be used to sling and swath an arm, to secure splints to an extremity or to perform in any other function requiring an elastic wrap.

Features:

- Suitable for applying immediate direct pressure to larger areas
- Oversized 12 in. x 12 in. sterile non-adherent pad
- Simple application that is easy to remember under combat/survival stress
- Omission of pressure bars and hooks to moderate the use of fine motor skills
- Vacuum-sealed, low-cube packaging with NAR®'s signature Red-Tip Technology® red tear notches

Specifications:

- Packaged Dimensions:
L 5.25 in. x W 6 in. x D 1.5 in.
- Bandage: L 70 in. x W 6 in.
- Pad: H 12 in. x W 12 in.
- Weight: 5 oz

Abdominal/Stump ETD®

| ITEM# | NSN# |
|---------|------------------|
| 30-0012 | 6510-01-541-8121 |



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**