

# EMERGENCY TRAUMA DRESSING - ETD®

8 X 10 IN. LARGE FLAT FOLDED ETD®

REV072823



## Rapid, controlled, effective bandaging with 8 x 10 in. pad ideal for larger wounds

Introducing the newest bandage in our proven series of Emergency Trauma Dressings (ETD®). This ETD® consists of a 4 in. resilient elastic wrap with sewn-in extra large 8 x 10 in. sterile non-adherent pad ideal for applying immediate direct pressure to larger areas. Engineered using data from over sixty-three after action reviews (AAR) of real-world battle dressing applications, the 8 x 10 in. ETD® also features a durable securing device and innovative Quick-Grip Roll Control configuration which facilitates a rapid, controlled, and extremely effective application process in the most demanding of conditions.

The ETD® is also suitable for use as a sling and swath for the arm, a secure fastener for an extremity splint or for any other function that requires an elastic wrap.

### Features:

- Suitable for applying immediate direct pressure to larger areas
- Extra large, sewn-in 8 x 10 in. sterile non-adherent pad
- Simple application that is easy to remember under combat/survival stress
- Omission of pressure bars and hooks to moderate the use of fine motor skills
- Vacuum-sealed, low-cube packaging with NAR®'s signature Red-Tip Technology® red tear notches

### Specifications:

- Packaged Dimensions: L 4.125 in. x W 7 in. x D 1 in.
- Bandage: L 70 in. x W 4 in.
- Pad: H 8 in. x W 10 in.
- Weight: 4 oz

### 8 x 10 in. Flat Folded ETD®

Item#	NSN#
30-0165	6510-01-654-6801



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**