

# EXPEDITIONARY CASUALTY RESPONSE KIT (ECRK)

REV081123

## Lightweight platform designed to rapidly access and organize all your essential trauma management equipment

The Expeditionary Casualty Response Kit (ECRK) features silent, dirt and sand-resistant zippers that close over a large main compartment with six (6) segregated removable equipment pouches; 4 Bleed Pouches, 1 Airway Pouch, and 1 Miscellaneous Pouch. Additionally, a separate lid compartment secures the R.E.E.L. Non-Rigid Litter while two internal pouches on the internal lid compartment and two external pouches on each side of the bag accommodate additional trauma management essentials.

### Features:

- Constructed of 500D Cordura® solution dyed nylon
- Sand- and water-resistant, self repairing zippers with silent 550 pull cords
- Separate lid compartment for the R.E.E.L. Non-Rigid Litter
- Six segregated removable pouches secured by hook and loop:
  - 4 Internal Bleed Pouches with Red pull handle
  - 1 Airway Pouch with Blue pull handle
  - 1 Accessories Pouch with Coyote pull handle
- 2 internal lid pouches to secure additional medical products
- 2 external side pouches on either side of the kit for additional medical products
- 2 external straps with dual adjust MILSPEC IR reduced buckles to secure the kit in a vehicle or boat

### Specifications:

- Case closed: H 17 in. x W 16.5 in. x D 10 in.
- Weight: 18 lb

OPTION	ITEM#	NSN#
Rx ONLY Kit	80-0948	----
Bag Only	80-1047	----

### Available Color:



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277



# EXPEDITIONARY CASUALTY RESPONSE KIT (ECRK)

REV081123

## ECRK Contents:

- 6 x Combat Application Tourniquet, (C-A-T®)
- 1 x Blizzard Blanket
- 4 x NAR® Pre-Lubricated Nasopharyngeal Airway, 28F
- 4 x ARS® Needle Decompression Device, 14g, 3.25 in.
- 4 x Combat Gauze®, Z-Fold Hemostatic Dressing
- 6 x Emergency Trauma Dressing, ETD®, 6 in.
- 10 x NAR® Compressed Gauze, 4.5 in. x 4.1 yd
- 4 x USMC HyFin® Chest Seal Combo Pk (1 HyFin® Vent & 1 HyFin® Chest Seal)
- 6 pr Bear Claw® Gloves, Large Tan
- 6 x Emergency Trauma Dressing, ETD® ABD 8 x 10 in.
- 1 x Minor BurnTec Kit
- 1 x NAR® Trauma Shears, Large, 7.25 in.
- 6 x Combat Casualty Response Cards
- 2 x Surgical Tape 1 in.
- 1 x R.E.E.L. Non-Rigid Litter
- 10 x Antiseptic Wipes
- 1 x Combat Casualty Reference Card
- 2 x SAM® Splint II
- 2 x Permanent Marker, Large Black
- 3 x Elastic Bandage 6 in.
- 2 x Triangular Bandage, 45 in. x 45 in. x 63 in.
- 2 x Mini Duct Tape, 2 in. x 100 in.
- 1 x Model 7 Rescue Hook
- 6 x Polycarbonate Eye Shield
- 1 x Expeditionary Casualty Response Bag



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**